SILVERGEAR

PULSE OXIMETER

HEART RATE & OXYGEN SATURATION

Manual

NL • Handleiding

FR • Guide

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IT • Manuale

ES • Manual

Thank you for purchasing our Pulse oximeter. This manual will provide instructions for effective and safe use of the product. Please read it carefully before use and store the manual in a safe place for future reference. If you ever pass the product on to others, this manual should be included.

1. Explanation of graphic symbols

Label	Explanation
	IMPORTANT Refer to instruction manual recommended.
	WARNING These warning notes must be observed to prevent any injury to the user.
	These notes give you useful additional information on the installation or operation.
	Indicates the item is a medical device.
	Degree of protection = IP22: Protected against solid objects over 12mm, e.g. fingers and against direct sprays of water up to 15 degrees from the vertical.
	Equipment classification: Type BF Applied Part
	Manufacturer
	Date of manufacture
	Serial number
	LOT number
	Authorized representative in the European community
	WEEE (Waste Electrical and Electronic Equipment) This product must not be disposed of together with domestic waste.
	CE symbol (conformity symbol)
	Unique Device Identifier

	Operating temperature
	Keep the device dry.

2. Safety instructions

- For intended use only, as described in this manual.
- This device is designed for general household use. Consult your doctor prior to use if you are concerned about health matters.
- Do not use to diagnose or treat any medical conditions.
- Do not use with medical implants such as pacemakers.
- Before use, inspect the device and its parts. Do not use in case of damage.
- Do not use the product in extreme hot, cold or humid environments.
- Do not place near heat sources such as radiators, cookers, etc.
- Do not place, store or use the product near water.
- Do not let moisture come into contact with any electrical parts.
- Do not attempt to repair, disassemble or alter the product yourself.

What is pulse oximetry?

Pulse oximetry is the most accurate non-invasive method for measuring blood oxygen levels. Normal blood oxygen levels in humans are considered between 95-100%. If the level is below 90%, it is considered low, resulting in hypoxemia with the blood not supplying enough oxygen to the body to function properly.

SpO2 classifications (oxygen saturation)

Note: The following table does not apply to people with certain pre-existing conditions such as asthma, heart failure or respiratory diseases, or at altitudes above 1,500 metres. If you have a pre-existing condition, always consult your doctor to evaluate your measurements.

SpO2 (oxygen saturation)	Classification - Action needed
99-95%	Normal - no action needed
94-91%	Decreased range - medical attention recommended
< 90%	Critical range - Urgent medical attention required


Note: Oxygen saturation levels decrease as altitude increases.

2. Package content

- Pulse oximeter
- Manual

3. Description of parts


1. Pulse bar
2. Power button
3. Low-battery indicator
4. Oxygen saturation (value in %)
5. Pulse rate (value in beats per minute)
6. Battery compartment

4. Operation

4.1 Installing batteries

1. Open the battery compartment on the bottom of the pulse oximeter.
2. Insert two 1.5V LR03 AAA batteries (not included), with regards to polarity (+ and -).
3. Close the battery compartment.

4.2 How to use


4.2.1 Preparations

Note:

- Wash your hands.
- Make sure your nails are free of nail polish or artificial nails.
- Rest for 15 minutes before measuring.
- Measurements should be taken in a quiet environment.
- Make sure you are relaxed and comfortable.
- Check that your hand is well circulated and warm to the touch.
- Sit relaxed but upright on a chair.
- Make sure your arm is supported and lower than your heart.
- Save the measured values for your GP just to be sure.

4.2.2 How to take a measurement

1. Carefully attach the pulse oximeter to the index finger.
2. Press the Power button to turn it on.
3. Sit still and relax your hand for about 10-15 seconds while taking a measurement.
4. Observe the values on the display. If necessary, short press the Power button to rotate the display (vertical/horizontal).
5. When you remove the device from your finger, it automatically switches off after about 15 seconds.

Note: Take three measurements every day spread throughout the day (morning, afternoon, evening) to get an accurate indication of your actual oxygen levels in your blood.

5. Maintenance

- Clean the product before and after each use with a cotton pad and some alcohol-based liquid (70%). Let the alcohol air dry to effectively kill germs.
- Remove the batteries before long-term storage.
- Store it in a dry, cool place away from direct sunlight.

6. Troubleshooting

Symptoms	Possible cause	Solution
No measurements are displayed	1. Finger is not correctly positioned between sensor and LED. 2. Blood oxygen levels are too low to be detected.	1. Attach the product correctly. 2. Make sure nothing is restricting the blood flow.
Values keep changing.	1. Finger is not correctly positioned between sensor and LED. 2. The user is moving.	1. Make sure the product is correctly attached. 2. Try to stay still during the measurement.
Unable to switch on.	1. Batteries are drained. 2. Batteries are installed incorrectly. 3. The product is defective or damaged.	1. Replace the batteries. 2. Remove and reinsert the batteries correctly. 3. Contact the distributor.
Display suddenly turns off.	1. The product switches off automatically after 15 seconds. 2. Batteries are drained.	1. This is normal. Just turn on the pulse oximeter again. 2. Replace the batteries.

7. Technical specifications

Battery	2x 1.5V LR03 AAA (not included)
Input	3V
Medical class	Ila
Measurement range	SpO2: 35-99% Pulse: 25-250 bpm
Accuracy	SpO2: $\pm 2\%$ for 70%-99% Pulse: ± 2 bpm
IP rating	IP22
Operating temperature	10~40°C
Relative humidity	15-95%RH
Atmospheric pressure	70~106kPa

8. Disposal

Dispose of the article and the packaging materials in accordance with current local regulations.

The crossed out wheellie bin symbol means that this product shall not be disposed of with normal household waste. Electronic and Electrical Equipment not included in the selective sorting process are potentially dangerous for the environment and human health due to the presence of hazardous substances. Please dispose of responsibly at an approved waste or recycling facility.

The CE marking is a certification mark that indicates conformity with health, safety and environmental protection standards for products sold within the European Union. The CE marking is the manufacturer's declaration that the product meets the requirements of the applicable EC directives.

Recycling helps to reduce the consumption of resources and reduce the environmental impact. If you wish to discard the packaging, comply to the environmental legislation in your country.

Karsten International B.V.

Overschiestraat 63, 1062 XD Amsterdam, The Netherlands
Tel: +31-(0)20-261 4883

Shenzhen AOJ Medical Technology Co.,Ltd.

Room 301 & 4F Block A, Building A, Jingfa Intelligent Manufacturing Park, Xiaweiyuan, Gushu Community, Xixiang Street, Bao'on District, 518126 Shenzhen, China.

Share Info GmbH

Heerdter Lohweg 83,40549 Düsseldorf, GERMANY
Tel: 0049 179 5666 508
E-mail: EU-Rep@share-info.com

Declaration


EMC of this product complies with IEC60601-1-2 standard.

The materials which the user can come into contact with have no toxicity and no action on tissues comply with ISO10993-1, ISO10993-5 and ISO10993-10.

Measurement Principle

Oxygenated blood absorbs light preferentially at 905nm (near infrared light), whereas deoxygenated blood absorbs light preferentially at 660nm (red light). A pulse oximeter works by passing a beam of red and infrared light through a pulsating capillary bed and then measure the amount of red and infrared light emerging from the tissues via a sensor. To improve accuracy, the Oximeter uses a proprietary algorithm to collect data from pulsatile arterial blood and excludes local noise from the tissues. The relative absorption of light by oxyhemoglobin (HbO) and deoxyhemoglobin is then calculated according to the Beer-Lambert's law and a quantitative measurement of the users' oxyhemoglobin status i.e. oxygen saturation level (SpO₂) is derived.

Due to the sensitivity of the pulse oximeter, finger should be kept stationary during measurement. It is recommended that you use this device for measurement before or after sports. Do not use for continuous monitoring.

Intended Purpose

The pulse oximeter is a reusable device and intended for spot-checking of oxygen saturation and pulse rate for use with the finger of adult patients.

Intended User

Pulse Oximeter could be used by Healthcare Professional and lay person.

Intended Patient

Pulse Oximeter is used to adult patients.

Indications

The pulse oximeter can effectively estimate and predict the blood oxygen concentration in the periphery of the body by non-invasive operation, so as to achieve the role of monitoring the changes of the respiratory disease with using the finger pulse oxygen as the monitoring index.

Contraindications

It is not intended to be used when the following happens:

pulse oximeter should not be fixed on the tissue injury site;

Not used for blood oxygen monitoring in hyperkinetic patients;

Not used for blood oxygen measurement under hypoperfusion.

Graphical Plots of Data Points


Figure 1 Bland-Altman graph for SpO2-SaO2


Figure 2 The regression line of the measured data

Assessing the Validity of a SpO2 Reading

You can check the quality of the pleth wave and the stability of the SpO2 values on the Patient Monitor to assess whether the sensor functions properly and whether the SpO2 readings are valid. Always use these two indications simultaneously to assess the validity of a SpO2 reading.

NOTE:

1. The SpO2 accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. Pulse oximeter measurements are statistically distributed, only about two-thirds of the measurements can be expected to fall within the specified accuracy compared to CO-oximeter measurements. The volunteer population in the studies composed of local healthy men and women from age 18 to 46, with variations of skin pigmentations. And the accuracy cannot be assessed by a function tester.
2. The pulse rate accuracy is obtained by comparison to the pulse rate generated with an arterial oxygen simulator (also an electronic pulse simulator).

3. Generally, the quality of the SpO2 pleth wave reflects the quality of the light signals obtained by the sensor. A wave of poor quality manifests a decline of the signal validity. On the other hand, the stability of the SpO2 values also reflects the signal quality. Different from varying SpO2 readings caused by physiological factors, unstable SpO2 readings are resulted from the sensor's receiving signals with interference. The problems mentioned above may be caused by patient movement, wrong sensor placement or sensor malfunction. To obtain valid SpO2 readings, try to limit patient movement.

Appendix 1 EMC Information

Guidance and manufacturer's declaration-Electromagnetic emission		
Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Model AOJ-70D should assure that it is used in such an environment.		
Emissions	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Pulse Oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	N.A.	N.A.
Voltage fluctuations/flicker emissions IEC61000-3-3	N.A.	

Guidance and manufacturer's declaration-Electromagnetic emission

Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Model AOJ-70D should assure that it is used in such an environment.

Emissions	IEC 60601 test level	Compliance level	Electromagnetic environment- guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable	
Surge IEC 61000-4-5	Not applicable	Not applicable	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable	
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC61000-4-6	Not applicable	Not applicable	
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	

NOTE: UT is the a.c. mains voltage prior to application of the test level

Guidance and manufacturer's declaration - electromagnetic Immunity

The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of Pulse Oximeter should assure that it is used in such an environment.

Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27	27
	450	430 – 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9	9
	745							
	780							
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28	28
	870							
	930							
	1720	1700–1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28	28
	1845							
	1970							

	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9	9
	5500							
	5785							

Guidance and manufacturer's declaration - electromagnetic Immunity				
Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IEC 60601-1-2 Test Level (A/m)	Compliance level (A/m)
	30 kHz	CW	8	8
	134,2 kHz	Pulse modulation 2.1 kHz	65	65
	13,56 kHz	Pulse modulation 50 kHz	7,5	7,5

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