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EN Users Manual HAPI-030D

1. Selecting custom-fit Ear Bud. The Hearing aid comes with 2 different size of ear buds. Before use, choose from either L or M size to attach to the ear piece as shown in Pic#1. It is important to choose the appropriate ear bud size for maximum comfort and sound quality.

Caution: Never insert the Ear Bud alone in your ear, it may be difficult to remove from the ear if it is not a attached aid.

2. Inserting Battery

- Open the battery door in a direction as shown in Pic#2 * Insert battery Pic#3
- Close the battery door

3. Wear the hearing aid in your ear by placing the Ear Bud in the ear & positioning the Ear Piece comfortably over the ear. To facilitate insertion of Ear Bud, lightly twist the unit left and right to achieve a snug fit. There is no need to push it into your ear strongly as shown in Pic#5

4. The hearing aid is triggered "ON" as soon as you closed the battery door, and volume is present level 4 as shown in Pic#4. To attain desired listening volume, 7 - level selection from 1 to 7 could be adjusted. Pic#4 Press "Tone" button to select among 4 Tones (Medium/High/Wide/Low) according to your personal hearing need as shown in Pic#6

5. Turn the unit "OFF" by when not in use by adjusting volume to minimum (level 0) or opening the battery door to the lock hinge as shown in Pic#7.

Cleaning your Hearing Aid. We recommend that you clean your hearing aid regularly (e.g. twice a week) to remove any earwax that may build up a result of regular use with a cotton tip or brush.

Caution: The hearing aid is not water proof product. Never immerse it completely in watery IEC 529 (ECC) level 0>. It is suggested to place the hearing aid back to the carrying case after daily usage and well store it at room temperature with good ventilation and humidity<30%, and avoid contact with corrosive gases.

Package Contents:



Hearing Aid
2 Ear buds (L/M)
3 Batteries (1,5V Alkaline Button Cell Battery LR44)
Carrying Case
Users manual

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