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PZN 11487206

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PZN 11487206

### ENGLISH

*Hosiery aid for compression hosiery (Instructions for use)*

#### Intended use

The donning aid is particularly useful for pulling on support and compression hosiery with restricted mobility.

#### Use

First set the hosiery aid down on a firm surface. The broad side faces you, whilst the open side of the semi-circular arc faces forwards.

Now pull the compression hosiery over the inner arc. The hosiery toe must be facing forwards here. Slide the upper rim of the compression hosiery over the arc edge until the toe of the compression hosiery appears at the upper edge of the arc.

Now place the hosiery aid on the floor and slide the foot, toes first, into the compression hosiery.

Slide the foot right down and check that the toe and heel sections of the hosiery are sitting correctly.

Once the hosiery is correctly positioned, hold the hosiery aid by the right and left frame handles and draw the hosiery aid in the direction of the knee. Upon reaching the knee draw the hosiery aid back and away to remove it.

You now need only pull up the top section of the hosiery. When doing so, ensure that the compression hosiery is not creased and that it cannot be pulled up any further.

#### Cleaning and care

The hosiery aid can be wiped with a damp cloth if necessary.

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**Dr. JUNGHANS Medical GmbH**