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
	FR	DE	NL	EN
①	SpO2 : saturation d'oxygène	SpO2: Sauerstoffsättigung	SpO2: zuurstofverzadiging	SpO2: oxygen saturation
②	Niveau de batterie	Batteriestand	Batterijniveau	Battery level
③	PRbpm : Rythme cardiaque	PRbpm: Herzrhythmus	PRbpm: hartritme	PRbpm: Heart rate
④	Bouton Marche/Arrêt	Ein/Aus-Taste	Aan/uit knop	On/Off button
⑤	Graphique du rythme cardiaque	Balkengrafik Pulsrhythmus	Staafigrafiek polsfrequentie	Heart rate graph
⑥	Intensité du pouls	Pulsstärke	Hartslag	Pulse intensity
⑦	PI	PI	PI	PI

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FR

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FINGERTIP PULSE OXIMETER

EN

Warning: Please read this instruction manual carefully before using the instrument and keep it for future use.

- Prolonged use of the oximeter may cause discomfort or even pain, especially in patients with circulatory problems. We do not recommend leaving the sensor on the same finger for more than 2 hours.
- More fragile patients should be monitored more carefully in the placement process. The device should not be placed over edema or sensitive tissue.
- The light (infrared is invisible) emitted by the device causes eye damage. Therefore, the user should not stare at this light under any circumstances. The patient must not have nail polish or any other type of cosmetic product on the fingers. The patient's fingernails should not be too long. Please read the sections on clinical restrictions and hazards carefully. This device is not a therapeutic device.

A. SAFETY INSTRUCTIONS

- This appliance is intended for home use only. Use the unit only as described in this manual.
- Never attempt to repair the unit yourself. Always have it serviced by a qualified service center or repair shop.
- Keep the unit away from heat sources such as radiators to avoid deformation of the plastic parts.
- Do not use abrasive cleaners to clean the unit.
- Do not use the pulse oximeter if it is damaged.
- This product was calibrated before leaving the factory.

Warning:

- Explosion Hazard - DO NOT use the oximeter in environments containing flammable gases (e.g., flammable anesthetic gases).
- DO NOT use the oximeter while the patient is undergoing MRI and CT scans.
- Allergy to rubber is a contraindication to the use of this device.

Warnings:

- Keep the pulse oximeter away from dust, vibration, corrosive substances, explosive substances, high temperatures, and moisture.
- Turn off the oximeter if it gets wet.
- Do not use the oximeter immediately after moving it from a cold to a warm or humid environment.
- DO NOT press the front panel buttons with sharp objects.
- Using the oximeter on a finger that is too thin or too cold may result in an inaccurate SpO2 and pulse measurement. Be sure to place the oximeter on a finger that is large enough, such as the thumb or middle finger, and insert it deep enough into the probe.
- Do not use the device on infants or newborns.
- The product can be used on children over 4 years of age and on adults
- The device may not work for all patients. If you are unable to obtain stable results, discontinue use.
- The data update time is less than 5 seconds and is variable depending on the pulse rate of the individual.
- If abnormalities appear on the screen during the testing process, remove your finger from the device and reinsert it to restore normal use.
- The normal useful life of the device is three years after it is first turned on.
- The instrument does not have a low battery alarm, only a low battery indicator. Please change the batteries when they are low.
- The oximeter does not have an alarm function for specific parameters. Do not use it in situations where an alarm is required.
- Remove the batteries if the unit is not used for more than one month. Otherwise, they may leak.
- A flexible circuit connects the two parts of the device. Do not twist the connection. Do not pull on it either.

B. HOW TO USE IT

Warnings:

1. If the finger is not inserted correctly, measurements may be inaccurate.
2. Excessive ambient light, such as from a fluorescent lamp, infrared heaters, direct sunlight, etc., may affect the measurement result.
3. Excessive patient movement or extreme electrosurgical interference may also affect measurement accuracy.
4. The patient should not have nail polish, false nails or any other type of cosmetic product on the fingers. This product is not suitable for continuous monitoring of patients. The results are distorted in case of carbon monoxide poisoning. Therefore, we do not recommend the use of this device in these circumstances.

Clinical Restrictions:

1. Since the measurement is taken on the basis of the arteriolar pulse, it is necessary that the subject's pulsatile blood flow power is sufficient. The SpO2 (plethysmography) waveform will be reduced in a subject with a weak pulse due to shock, low room/body temperature, severe bleeding, or use of vasoconstrictors. In this case, the measurement will be more sensitive to interference.
2. The presence of intravascular dyes in some patients (methylene blue, indocyanine green, and acid indigo blue), carboxyhemoglobin (COHb), methemoglobin (MetHb), or sulfhemoglobin, or the presence of jaundice may distort SpO2 determinations.
3. Medications such as dopamine, procaine, prilocaine, lidocaine, and butacaine may also be a factor in causing significant errors in SpO2 measurements.
4. Because the SpO2 value serves as a reference value for the assessment of anemic anoxia and toxic anoxia, some severely anemic patients may have good SpO2 measurements.

1) Insert and/or replace the batteries

- Open the battery cover on the back of the unit.
- Insert two "AAA" batteries. Please observe the polarity of the batteries.
- Replace the battery cover.

2) Operation

- Open the oximeter clamp by pressing firmly on the side opposite the ON button and insert your index finger or thumb into the appropriate position. Release the clamp.
- Make sure your finger is fully inserted. (When your finger is inserted into the unit, the surface of your fingernails should be facing the same side as the display).
- Press the On/Off button once ④.
- Stay still and avoid applying pressure to the device tightly during use.

O 8690 / xxx

- Wait a few seconds for the measurements to appear on the display.

Do not shake the finger and make sure the patient is comfortable during the process. Body movement is not recommended during the measurement.

Look at the screen to get the information.

If necessary, press the ON button briefly to select one of the 6 available display modes (rotating the screen in case a person accompanying you wants to be able to read the displayed values for example).

The ① SpO2 value corresponds to the level of oxygen contained in the red blood cells after they have passed through the lungs. It is therefore primarily to verify that a sufficient amount of oxygen in the red blood cells. Its result is expressed on the oximeter in percentage (between 95 and 100% for a normal rate). When it is less than 90% the person is said to be in desaturation.

The value ③ PRbpm corresponds to the pulse rate otherwise known as heart rate, the PR calculates the number of heartbeats per minute. This is an important index of your physical condition and is expressed as a whole number. Between 50 and 90 beats per minute, we can say that the heart rate is normal and common to most people.

The ⑦ value is the perfusion index (PI) or pulsatile signal intensity: Indicates the signal intensity of the arterial pulse. As a curve or as a vertical bar, higher or lower depending on the signal strength.

- When you are finished, remove your finger from the oximeter. The oximeter will automatically turn off after about 10 seconds.

3) Changing the settings

This product was calibrated before leaving the factory.

It is not necessarily necessary to modify the settings. However, you can, at each use, choose to change the settings (alarm function in case of measurement below or above a defined threshold the value will flash on the screen to signify the alarm)

Press the On/Off button ④ once to turn on the device and wait for it to turn on. Then press for at least 5 seconds to enter the device settings.

Press the ④ button once to select the line of the setting you want to change. When the line you want to change is in yellow, long press to start the adjustment: the values will increase slowly. Release the button when the value you wanted to select is reached.

Repeat for all the setting values you want to change.

When you are done, go to "exit" and long press the ④ button to exit.

C. CARE AND MAINTENANCE

Do not use abrasive cleaners to clean the unit.

Do not immerse the oximeter in any liquid.

Wipe the surface with a soft cloth dampened with medical alcohol.

Do not spray liquids directly onto the unit.

If you use water to clean the oximeter, the temperature should be below 60°C.

Do not use disinfection

D. TROUBLESHOOTING

Problem	Possible cause	Solution
Unable to display spo2 and pulse correctly	<ol style="list-style-type: none"> 1. The unit will automatically shut down if it does not receive a signal for 5 seconds. 2. The patient's SpO2 is too low to be detected 	<ol style="list-style-type: none"> 1. Place the finger correctly and repeat the procedure. 2. Try again; go to the hospital for diagnosis if you are sure the device is working properly.
the display of the spo2 and pulse is not stable	<ol style="list-style-type: none"> 1. The finger is not inserted to a sufficient depth. 2. The finger or the patient moves. 	Place the finger correctly and repeat the procedure. Wait for the patient to come to rest.
impossible to turn on the device on	<ol style="list-style-type: none"> 1. The batteries are dead or very weak. 2. The batteries are not inserted correctly. 3. Device malfunction 	<ol style="list-style-type: none"> 1. Change the batteries. 2. Reinstall the batteries. 3. Please contact your local service center.
The screen goes out suddenly	<ol style="list-style-type: none"> 1. The unit will automatically shut down if it does not receive a signal for 5 seconds. 2. The batteries are nearly exhausted. 	<ol style="list-style-type: none"> 1. Normal. 2. Change the batteries.

E. FEATURES

ABS and silicone materials

Power supply: 2 AAA batteries not included

SpO2 measurement range: 70% - 99

Pulse measurement range: 30 bpm - 240 bpm

Measurement accuracy: $\pm 2\%$ for SpO2 values from 70 % to 100 %, and not specified for values below 70 % . ± 2 bpm or $\pm 2\%$ (whichever is greater) for pulse

Measurement performance under low irrigation conditions: SpO2 and pulse can be displayed correctly when the pulse-irrigation ratio is 0.4%. The error for SpO2 is $\pm 4\%$, the error for pulse is ± 2 bpm or $\pm 2\%$ (whichever is greater).

Resistance to ambient light: the difference between the value measured under artificial light or natural light in a closed environment and in a dark room is less than $\pm 1\%$.

The unit shuts down after 8 seconds of inactivity.

The European directive 2012/19/EC about electronic and electric waste, requires that you can't throw away defective domestic appliance with common waste. Used device has to be collect separately to optimize the recuperation rate and the recycling of the materials to reduce the impact on the health and the environment.

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capacities, or by persons lacking in experience or knowledge, unless they are supervised or instructed on the use of this apparatus by a person responsible for their safety. Children should be supervised to ensure they do not play with this appliance.

Precision and precautions about the battery use

- Throwing away a used battery in the environment or with household waste pollutes and avoids the recuperation of reusable materials.
- It is important to limit the battery consumption and to respect the following rules:
 - Favor alkaline battery (last longer than saline battery) and when it's possible use rechargeable battery.
 - Leave used batteries and accumulators in adapted containers that you can find in some retailers. In that way, metals don't pollute the environment and the health.
- Battery has to be inserted by respecting the polarity show on the apparel and the battery. An incorrect positioning can damage the apparel, cause leaking and even provoke fire and the explosion of the battery.
- To assure a good functioning, battery has to be in a good condition. If there is a problem, put new battery.
- Never try to recharge battery not made for this. It could leak, warm up, cause fire or explode.
- Replace all the battery at the same time. Never mix up alkaline battery with saline battery.
- Used battery has to be removed.
- Also, remove the battery if you're not using the apparel for a long time.
- Never try to short-circuit the terminal of a battery.
- Never throw away the battery in fire because they can explode.
- Take out the accumulator from the apparel before reload it.
- We recommend that an adult watch if a child change the battery to respect those rules.
- If you swallow a battery, go immediately to your doctor or the hospital.
- For the best ongoing performance, we recommend using only fresh, high-quality, alkaline batteries.

Importé par / Importiert von / Geïmporteerd door / Imported by :
CAPTELEC, 59170 Croix - FRANCE