Relax Vital

Beintrainer

Deutsch

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Relax Vital Beintrainer

English

Instruction Manual
from page 6

Relax Vital leg trainer

Français

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Congratulations on your purchase of Relax Vital leg trainer.

Millions of people live a sedentary lifestyle. Immobility contributes to muscle atrophy. Many careers require the individual to remain seated for many hours per day. This immobility can contribute to multiple health concerns.

Relax Vital leg trainer is a revolutionary alternative that aids in problems caused by lack of mobility and/or an inactive lifestyle.

Relax Vital leg trainer can improve our blood circulation.

Relax Vital leg trainer is a passive gym that we can take advantage of while seated at work during long meetings, watching TV, reading, etc. without any physical effort on your part, thanks to its motorized continuous reciprocating movement.

Relax Vital leg trainer is a revolutionary exercise machine that moves your legs back and forth, alternatively, and smoothly massages your feet while you remain seated for long periods of time. This stimulates blood circulation and produces a sensation of comfort. It aids in back and leg pain, it could help decrease the proliferation of varicose veins, cellulite and could help prevent stiff knees.

Relax Vital leg trainer was designed with Ergo-Therapy in mind where its movements are precise, continuous and without vibrations. This makes the treatment very comfortable and smooth to the ankle and knee joints.

If you have any questions, please contact customer service via our e-mail address:

kundenservice@hsphanseshopping.de

HOW DOES IT WORK?

When we spend most of our day sitting in a chair, in your office, meetings, watching TV at home, etc., without moving the legs, poor blood circulation can occur. Blood accumulates in the legs, rather than circulating regularly. Relax Vital leg trainer works by moving the legs back and forth alternatively without any effort on your side so that blood circulates normally every time.

WHEN WOULD I SEE RESULTS?

Results are almost immediate. Relax Vital leg trainer maintains your legs moving while assisting your blood circulation, this in turn will prevent your legs from feeling tired.

HANDLING

Place the machine flat on the floor in front of the chair you are going to sit on. For your comfort, you can choose to leave it flat, or pull out the rear supports depending on your seating arrangement. Plug it into the nearest electrical outlet. Make sure the machine is in a cool and ventilated area.

For internal use only.

INSTRUCTIONS OF USE

1. Unit's location: must be flat on the floor. It has two rear legs. Pull out if necessary. **We recommend using some level of inclination, it helps with posture.**
2. Plug in the adapter to the nearest electrical outlet.
3. Plug in the other end of the adapter to the unit in the front of the machine.
4. Select the desired speed:
 1. **Normal Speed:**
recommended for continuous and normal use.
 2. **Fast Speed:**
when you plan to use the machine for short periods of time or less than an hour.
5. Relax Vital leg trainer has being designed to be used while seated only. **Never stand up on the machine**, as severe personal injury may result, and/or the weight of your body may damage the machine's internal mechanism and its functions.
6. Only use the machine while seated. Put your feet over both platforms. Do it in a relaxed way without doing any type of effort.
7. Enjoy while you work, read, watch TV, have lunch or dinner and many more.

SUGGESTIONS FOR BETTER RESULTS

1. Use in a cool and ventilated area.
2. Unplug the unit while not in use.

Due to our company's continuous improvement program, the images, videos or photos may vary from the final product.

CLEANING AND MAINTENANCE

Relax Vital leg trainer does not need any type of maintenance and it is designed for long periods of use. To clean it, turn the unit off and unplug it. Then wipe clean with a dry cloth.

WARNINGS

- This unit should not be exposed to water or splashes of water. Objects that contain liquids like a flower vase, coffee cup, etc., should not be near the unit.
- If you use an extension, keep the strip away from humidity or damaged cables.
- Keep fingers away from moving parts at all times.
- **Keep out of reach of children.**

DISPOSAL

Dispose of packaging material in an environmentally friendly manner and take it to the recycling collection.

This device is subject to the European Directive 2012/19/EU for Waste Electrical and Electronic Equipment (WEEE). Do not dispose of this device as normal household waste, but in an environmentally friendly manner via an officially approved disposal company.

Accumulators and batteries must be removed and disposed of separately from the device before disposal. In the interests of environmental protection, batteries and accumulators must not be disposed of with normal household waste, but must be handed in at appropriate collection points.

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HSP Hanse Shopping GmbH