Figure 1
Gráfico 1
Graphique 1
Bild 1
Imagine 1
Gráfico 1
Afbelding 1

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VELFORM® SATIN is a registered EU/CTM trade mark. International patents pending. All imitations will be prosecuted.

ENGLISH

Dear Customer:

Thank you for purchasing the Velform® Satin foot care system. This device runs on batteries (not included) and cares for your feet without damaging the skin. We hope you get the most out of your Velform® Satin. The device is recommended for one person only. For hygiene reasons, the roller should not be shared with others.

General Description

- A: Roller head
- B: Release button
- C: Safety lock
- D: Sliding ON/OFF switch
- E: Battery compartment cover
- F: Battery cover lock
- G: Roller for sensitive skin

Important:

1. This device is designed for use only on corns and calluses. Do not use it on uncallused skin, as this may cause injury.
2. Use the device only for purposes described in this instruction manual.
3. Each roller is recommended for one person only. For hygiene reasons, the roller should not be shared with others.
4. Do not use the device if you have diabetes or poor circulation.
5. If you have a skin condition or other medical problem, consult your physician before using the device.
6. Do not use on the skin for more than 2 to 3 seconds at a time.
7. Do not use on skin that is red, irritated, inflamed, infected or blistered.
8. Do not use on warts or open sores.
9. Discontinue use if you experience any pain or irritation.
10. Overusing the device on a single spot may cause skin irritation.
11. Seek medical attention if the skin becomes severely irritated or bleeds.
12. This device is designed for household use only.
13. This product was not designed to be used by anyone (including children) with reduced physical, sensory or mental abilities, or by people without adequate experience and/or knowledge, unless they are supervised by someone responsible for their safety or are given instructions on how to use the device.

Warnings

14. Do not place or store the device in direct sunlight. Store the device in a cool, dry place at a temperature between 0°C and 40°C. Keep the device dry when it is not in use.
15. Do not turn on the device if the roller is not in place.
16. Do not drop or insert anything into any opening on the unit.
17. Do not use the device outdoors.
18. To prevent accidents and injury to yourself or damage to the unit, while the device is in use, keep it away from the hair on your scalp and your eyebrows and eyelashes, as well as your clothing, and any wires, cables, shoelaces, etc.
19. Do not continually rub the unit on the same spot or apply strong pressure.
20. Do not disassemble the device.

Battery Installation

1. Remove the battery cover (battery cover lock) by pressing the tab.
2. Insert 2 AA alkaline batteries in the battery compartment following the battery direction sign printed inside the compartment.
3. Close the battery cover until it clicks into place and locks.

INSTRUCTIONS FOR USE

This device is designed to smooth rough skin on the feet and heels. Do not use it on any other part of the body.
 1. Wash the skin and make sure it is completely dry. Warning: Do not use the device on broken or bleeding skin.
 2. Make sure the roller head is in place on the unit. See the Changing the Roller Head section on how to remove and replace the roller.
 3. Turn on the device by pressing the safety lock in the centre of the ON/OFF switch and sliding the button up.
 4. Gently slide the device (back and forth or side to side) over the callus/corn for 2 to 3 seconds. Do not press down hard. Gently glide the device over the area to be treated. The unit will stop if it is pressed down too hard on the skin. (See Figure 1)
 Stop and check if you have achieved the desired softness. If not, go over the area for 2 or 3 more seconds and check again. Continue using this method until you have achieved the desired softness.
 5. Turn off the device. If the safety lock is on, unlock it by pressing the centre of the ON/OFF button and then sliding the button down.
 6. Rinse the skin or use a damp towel to remove the dead skin particles.
 7. To clean the device, consult the Cleaning section.

Changing the Roller Head

The roller head should be changed when necessary for best results. The roller head will last for several sessions, depending on the frequency of use and the size of area being treated. Replace it as needed. Make sure the device is turned off before replacing the roller head.

1. Hold the roller head with one hand and use the other to press the release button on the side of the device to release the roller head. The roller head will pop out. (See Figure 2)
2. Place a new roller head in the device and press down until you hear it click into place. Ensure the new roller head is firmly installed in the device. Warning about changing the roller head: Do not turn on the device if the roller head is not in place.

Cleaning

Warning: For best performance, clean the roller head after each use. The roller is water resistant so it can be rinsed in water.

1. Remove the roller head from the device and rinse it under running tap water.
2. Dry it with a clean, dry cloth.
3. Insert the roller head back into the device.

Storage

When the device is not in use, keep it in a dry place out of reach of children.

Do not store the device in an area with excessive heat or moisture.

CERTIFICATES AND DIRECTIVES

Velform® Satin complies with the following EC directives:

- Electromagnetic Compatibility Directive (EMC) 2004/108/EC
- RoHS Directive 2011/65/EC

Warranty

This product is covered by a warranty against manufacturing defects subject to the time periods stipulated by the legislation in force in each country.

Technical Data: 2 1.5V AA batteries (not included)

Level of protection against water: IPX2

CAUTION

DISPOSAL OF MATERIALS

 A symbol of a crossed-out wheeled bin means you should find out about and follow local regulations about disposing of this kind of product.

Do not dispose of this product as you would other household waste. Dispose of this device in accordance with the corresponding local regulations.

Electrical and electronic devices contain hazardous substances that can have harmful effects on the environment and/or human health and should be recycled properly.

Australia & New Zealand only

This product has the benefit of certain Consumer Guarantees. These are prescribed by the Australian Consumer Law & the New Zealand Consumer Guarantees Act 1993 both of which provide protection for consumers. There is no express warranty for this product in Australia or New Zealand. The above paragraph refers to other countries.

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