Mannsberger®

Dual Bikeý
Instructions

Dual Bikeý
Handleiding

Dual Bikeý
Instructions

Dual Bikeý
Gebrauchsanweisung

Istruzioni

Instrucciones

Precautions

- Read instructions thoroughly before first use and keep them safe for future reference.
- Only operate the Dual Bike on level ground. Inspect the exercise bike for worn or loose components before each use. The safety and integrity of the Dual Bike can only be maintained when the unit is regularly inspected for damage. It is the sole responsibility of the operator/owner to ensure that regular maintenance is performed. Do not operate unit with worn or damaged components. Use only manufacturer supplied components to maintain/repair unit.
- Only operate the unit in a sitting position on a sturdy chair. DO NOT operate the exercise bike in a standing position. Consult your physician before beginning any exercise program. If you experience dizziness, faint, or ill during exercise, stop exercise immediately and consult your physician.
- Do not wear loose clothing while operating the exercise bike.
- Use caution when mounting and dismounting the exercise bike.
- This exercise bike is for consumer and home use only.
- The Dual Bike is not a medical device.
- Maximum user weight: 220 lbs. (100 kg.)
- Store unit in a cool and dry place



Warning

- Keep exercise bike away from children.
- Avoid contact with spinning pedals as they can cause injury.


- Gradually reduce pedaling speed when coming to a stop.
- The exercise bike should only be used after a thorough review of the instructions.

Parts List

- | | |
|---|-----------------------------|
| A x 4: Small screws & washers | L: Main frame tension knob |
| B x 2: Front side base spacers | M: Digital Exercise Counter |
| C x 2: Back side base spacers | N: Snap buttons |
| D x 1: Screw size 3.75" with washer and nut | O: Magnetic sensor |
| E x 1: Screw size 3.50" with washer and nut | P: Front base tabs |
| F x 1: wrench tool | Q: Back base tabs |
| G x 1: Screwdriver hex key tool | |


Use screw size 3.75"(D)


Use screw size 3.5"(E)

1. Insert the front base spacers (B) between the front base tabs (P) of the main frame (K). Then secure with screw size 3.75", washer, and nut (D). Tighten screw with hex key (G) and wrench (F) included.
2. Repeat step 1 process for the back base (Q) and frame (K), using the back side base spacers (C), screw size 3.5", washer, and nut (E) as shown.

3. Insert the pedals (I) on the underside of main frame (K) as shown. Secure with small screws and washers (A) using hex key (G). Note: make sure the magnet (J) is facing towards the magnetic sensor (O) on the main frame (K) as shown. Increase or decrease tension by turning tension knob on the pedals (I).

4. Insert handlebar (H) into the top main frame (K). Secure with small screws and washers (A) using hex key (G). Increase or decrease tension by turning tension knob on the handlebar (H).

5. Turn main frame knob (L) and push the snap buttons (N) to adjust the desired angle and height. Then tighten knob to secure.


How to Operate

6. Sit on a sturdy chair. Insert your feet into the pedals and grab the handlebars and rotate your arms and legs as desired. The digital counter will automatically turn ON and turn OFF when you start or stop rotating the pedals or if the MODE key (M) is not pressed within 4 minutes.

7. Digital Exercise Counter Specifications:

- | | |
|--------|--|
| CNT | = counts number of pedal resolutions |
| 0-9999 | |
| TM | = counts the time 0-99:59 |
| CAL | = counts the calories 0-99.99 |
| DIST | = counts the distance 0-99.99 |
| TCNT | = counts total of pedal resolutions 0-9999 |
| SCAN | = automatically scans through all counts every 3 seconds |

- Press the MODE red button (M) at any time to select any counter specification above.
- To reset the counter, press and hold the MODE key (M) for 3 seconds.

8. To replace battery, gently pull on the bottom of the digital counter face plate and slide the battery upward. Replace with LR44H button cell battery with the polarity (+) facing upward. Replace the digital counter face plate.

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Mesa Products International B.V.
Havenstraat 57 a, 7005 AG Doetinchem – NL